

Product Name:	Closed Liquid Transfer System of Square Storage Bottles
Product Application:	NEST launches a new version of liquid transfer system to provide an aseptic transfer solution for pharmaceutical, biotech and laboratories, enabling closed injection of liquid. Standard Luer head is used in the joint, making the operation easier and more convenient.
Product Materials:	PETG(bottle) + PP(cap) + TPE thermoplastic tube Meets <i>USP, Class VI standards</i>
product Specifications:	Temperature range: stored at room temperature Shelf life: 3 Years after date of production (ensure package is in good)
Sterilization:	Yes. Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
Pyrogens:	Non-Pyrogenic
RNase/DNase:	DNase/RNase free
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Sealing test: qualified Assembly inspection: qualified Freeze testing: qualified Capacity test Packaging inspection: qualified
Features	<ol style="list-style-type: none"> 1. TPE material of tubes is in line with the requirements for industrial pharmaceutical enterprises, and has passed the dissolution test carried out by a certified third party according to the complete BPOG dissolution test criteria. 2. Standard Luer head is used in the joint, making the operation easier and more convenient. 3. The transfer cap comes out of integrated injection molding, whose congruence guarantees its adaptation to other brands. 4. PETG material of the wall is of high transparency, strong mechanical strength, resistant to low temperature and convenient for observation and transportation 5. The unique thick collar design of bottle neck makes inside operation easier. 6. DNase/RNase/Pyrogen free. 7. Complete varieties comprise full specifications.

Product Range

PETG Square Storage Bottle with Tube

Cat. No.	Product Description	/Pack	/Case
353551 C50921-BHB060A	250mL PETG Square Storage Bottle with Tube, TPE Tube (60cm, 1/4" ID, 3/8" OD), Female Luer Lock Connector with Luer Plug, Vent Filter (0.22 µm, Φ33mm)	1	10
333551 C51021-BHB060A	500mL PETG Square Storage Bottle with Tube, TPE Tube (60cm, 1/4" ID, 3/8" OD), Female Luer Lock Connector with Luer Plug, Vent Filter (0.22 µm, Φ33mm)	1	10
334551 C51121-BHB060A	1000mL PETG Square Storage Bottle with Tube, TPE Tube (60cm, 1/4" ID, 3/8" OD), Female Luer Lock Connector with Luer Plug, Vent Filter (0.22 µm, Φ33mm)	1	10
352552 C50921-BBB060A	250mL PETG Square Storage Bottle with Tube, TPE Tube (60cm, 1/4" ID, 3/8" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 µm, Φ33mm)	1	10
333552 C51021-BBB060A	500mL PETG Square Storage Bottle with Tube, TPE Tube (60cm, 1/4" ID, 3/8" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 µm, Φ33mm)	1	10
334552 C51121-BBB060A	1000mL PETG Square Storage Bottle with Tube, TPE Tube (60cm, 1/4" ID, 3/8" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 µm, Φ33mm)	1	10
352553 C50922-AGB060A	250mL PETG Square Storage Bottle with Tube, TPE Tube (60cm, 1/8" ID, 1/4" OD), Female Luer Lock Connector with Luer Plug, Vent Filter (0.22 µm, Φ33mm)	1	10
333553 C51022-AGB060A	500mL PETG Square Storage Bottle with Tube, TPE Tube (60cm, 1/8" ID, 1/4" OD), Female Luer Lock Connector with Luer Plug, Vent Filter (0.22 µm, Φ33mm)	1	10
334553 C51122-AGB060A	1000mL PETG Square Storage Bottle with Tube, TPE Tube (60cm, 1/8" ID, 1/4" OD), Female Luer Lock Connector with Luer Plug, Vent Filter (0.22 µm, Φ33mm)	1	10
352554 C50922-AAB060A	250mL PETG Square Storage Bottle with Tube, 60 cm TPE Tube (60cm, 1/8" ID, 1/4" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 µm, Φ33mm)	1	10

Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
 Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
 Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
 NEST scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

ISO9001 & ISO13485 & ISO11137

333554 C51022-AAB060A	500mL PETG Square Storage Bottle with Tube, TPE Tube (60cm, 1/8" ID, 1/4" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 µm, Φ33mm)	1	10
334554 C51122-AAB060A	1000mL PETG Square Storage Bottle with Tube, TPE Tube (60cm, 1/8" ID, 1/4" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 µm, Φ33mm)	1	10

PETG Square Storage Bottle Transfer Cap

Cat. No.	Product Description	/Pack	/Case
353951 C50921-BHB060B	250mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60cm, 1/4" ID, 3/8" OD), Female Luer Lock Connector with Luer Plug, Vent Filter (0.22 µm, Φ33mm)	1	10
333951 C51021-BHB060B	500mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60cm, 1/4" ID, 3/8" OD), Female Luer Lock Connector with Luer plug, Vent Filter (0.22 µm, Φ33mm)	1	10
334951 C51121-BHB060B	1000mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60cm, 1/4" ID, 3/8" OD), Female Luer Lock Connector with Luer plug, Vent Filter (0.22 µm, Φ33mm)	1	10
352952 C50921-BBB060B	250mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60 cm, 1/4" ID, 3/8" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 µm, Φ33mm)	1	10
333952 C51021-BBB060B	500mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60 cm, 1/4" ID, 3/8" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 µm, Φ33mm)	1	10
334952 C51121-BBB060B	1000mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60cm, 1/4" ID, 3/8" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 µm, Φ33mm)	1	10
352953 C50922-AGB060B	250mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60cm, 1/8" ID, 1/4" OD), Female Luer Lock Connector with Luer Plug, Vent Filter (0.22 µm, Φ33mm)	1	10



Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
 Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
 Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
 NEST scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

ISO9001 & ISO13485 & ISO11137

333953 C51022-AGB060B	500mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60cm, 1/8" ID, 1/4" OD), Female Luer Lock Connector with Luer Plug, Vent Filter (0.22 μm, Φ33mm)	1	10
334953 C51122-AGB060B	1000mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60cm, 1/8" ID, 1/4" OD), Female Luer Lock Connector with Luer Plug, Vent Filter (0.22 μm, Φ33mm)	1	10
352954 C50922-AAB060B	250mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60cm, 1/8" ID, 1/4" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 μm, Φ33mm)	1	10
333954 C51022-AAB060B	500mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60cm, 1/8" ID, 1/4" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 μm, Φ33mm)	1	10
334954 C51122-AAB060B	1000mL PETG Square Storage Bottle Transfer Cap, TPE Tube (60cm, 1/8" ID, 1/4" OD), Male MPC connector with Female Sealing Cap, Vent Filter (0.22 μm, Φ33mm)	1	10

Technical Drawing



Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
 Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
 Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
 NEST scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)