

NEST[®]

Biofactory™
(Cap with Overcap)




+86-510-6800 6788

info@nest-wuxi.com

www.cell-nest.com

Biofactory™ (Cap with Overcap)



● Features

Vented Overcap and Solid Overcap can be used directly on the adaptor caps to block dust and bacteria. Vented Overcap is equipped with a 0.22 µm hydrophobic breathable membrane to achieve aseptic ventilation during liquid transfer.

Cat.No.	Description	/Pack	/Case
771322	NEST Biofactor 10 Chamber with 2 Solid Overcaps, TC, Sterile (with 12 Sterile Vented Overcaps Included Separately)	1	6
771422	Nest Biofactor 40 Chamber with 2 Solid Overcaps, Total Culture Area: 25295 cm ² , TC, Sterile, 1/pk, 2/case, (with 4 Sterile Vented Overcaps Included Separately)	1	2



740913



740901



740213

Cat.No.	Description	/Pack	/Case
740913	Vented Overcap for BioFactories, Adaptable for 740201, Individually Wrapped, Sterile	1	20
740901	Solid Overcap for BioFactories, Adaptable for 740201, Individually Wrapped, Sterile	1	20
740213	Adaptor Cap(740201) with Vented Overcap for BioFactories, Individually Wrapped, Sterile	1	10