Easy injection in three steps



Disposable **Pen Injector**



NEST Disposable Pen Injector is a type of contains a liquid medicine solution within a syringe. The medication is conveniently placed in a cartridge within the pen, eliminating the need to extract liquid medication, simply remove the pen cap, adjust the dial to the desired dose, and then

drugs, etc.










(♪) +86-510-6800 6788


(☑) info@nest-wuxi.com

Satisfy customer's needs



Product Feature

- Highly customizable: appearance, color, dosage, scale lines customizable.
- Currently available in 36/ 37 /50 /60 /74 /75 /80 units, which can be customized for different indications.
- Pressing force: 36/37 units <5N, 74/80 units <8N.
- Accurate dosage, professional R&D team and data validation support.
 Completed the FTO review in China and the United States, no patent
- risk.
- Applicable to liraglutide, simethicone, follicle stimulating hormone (FSH) drugs, insulin degludec and other drugs.
- Successfully applied for CDE number; CE MDR, FDA (510k) has been submitted for review.



Simple operation

It is easy for beginners to use, eliminating the need for professional guidance.







- 2、Dose regulation
- 3、Injection

Disposable Pen Injector Customization Solution

Product Name	Max. Single Dose	For Cartridge	Min. Single Dose	Usual Dose	Usual Scale Units	Applicative Drug
Disposable Pen Injector	36U	3mL	0.0208mL	0-12.5-25450	IU	Follicle stimulating hormone,FSH
	37U	3mL	0.01mL	0-0.25-0.5	mg	Semaglutide, 1.34mg/mL , 0.5mg
	50U	3mL	0.01mL	0-0.6-1.2-1.8-2.4-3.0	mg	Liraglutide
	60U	3mL	0.01mL	0-260	IU	Insulin Aspart series; Insulin glargine series
	74U	3mL	0.01mL	0-0.25-0.5-1.0 ;Customized	mg	Semaglutide(0.68mg/mL、1.34mg/mL、2.68mg/ml)
	75U	3mL	0.01mL	0-1mg; 0-1.7mg; 0-2.4	mg	Semaglutide
	80U	3mL	0.01mL	0-280	IU	Insulin degludec/insulin aspart,IDegAsp
SP style Disposable Pen Injector	37U	3mL	0.01mL	0-0.25-0.5	mg	Semaglutide, 1.34mg/mL , 0.5mg
	50U	3mL	0.01mL	0-0.6-1.2-1.8-2.4-3.0	mg	Liraglutide
	60U	3mL	0.01mL	0-260	IU	Insulin Aspart series; Insulin glargine series
	72U	3mL	0.01mL	0-450 IU	IU	Follicle stimulating hormone,FSH
	74U	3mL	0.01mL	0-0.25-0.5-1.0 ;Customized	mg	Semaglutide(0.68mg/mL、1.34mg/mL、2.68mg/ml)
	75U	3mL	0.01mL	0-1mg; 0-1.7mg; 0-2.4	mg	Semaglutide
	80U	3mL	0.01mL	0-280	IU	Insulin degludec/insulin aspart,IDegAsp
Disposable fixed dose Pen Injector	8U	1.5mL	0.08mL	0- inject	IU	Parathyroid Hormone ; Drug for Osteoporosis
	8U	3mL	0.08mL	0- inject	IU	Drug for Osteoporosis:Teriparatide